


VIVASALT *Mediteran*



INDEX

1. Composition of VIVASALT Mediteran.....	3
2. How VIVASALT Mediteran works?.....	3
3. Installation VIVAsalt Mediteran	4
3.1 Installation environment	4
3.2 Installation.....	4
4. Working with VIVASALT Mediteran.....	6
4.1 Salt.....	6
5. Cleaning and maintenance of VIVASALT Mediteran	9
5.1 Cleaning procedure	9
6. Solving Problems	13
7. Warranty terms	14
7.1 Contraindications (responsibility to patients).....	15
8. Storage of the device.....	15





Dear customer.

Congratulations on the purchase of a technologically advanced device for salt therapies of the VIVASALT brand. The device can serve you for a very long time if it is regularly maintained and serviced in accordance with these operating instructions, so we recommend that you always have it at your fingertips.

The instructions cover all information regarding installation, maintenance, troubleshooting, enforcement of the warranty and accompanying documentation.




1. Composition of VIVASALT Mediteran

- housing with electronics
- touchscreen
- fan
- mill
- power cable
- deflector.

2. How VIVASALT Mediteran works?

The VIVASALT Mediteran device is used only for the production of a dry aerosol of a medical salt that is blown in through the opening in the wall to the room. Typically, the device is located on the outer wall of the room into which the aerosol salt is blown in. This is called halotherapy or salt therapy.

A medical salt is put into the mill (the maximum particle size of the salt is 1 mm). The mill grinds the medical salt into a micron-size aerosol, which is blown with the fan into the room. Only aerosol is blown into the room, the larger pieces of salt hit the deflektor and fall to the ground.



3. Installation VIVAsalt Mediteran

Only an authorized person can mount and connect the device. Before installing, make sure that the person who installs the device will be authorized by the manufacturer VIVA SALT d.o.o., because in the event of a malfunction due to incorrect assembly or incorrect electrical connection, we do not accept warranty.

3.1 Installation environment

The environment in which VIVASALT Mediteran is mounted must be dry and warm.

Border range:

- the temperature between 18 and 25 degrees Celsius
- Moisture between 20% and 50%

TECHNICAL INFORMATION:

VOLTAGE: 230 V

WEIGHT: 7 Kg

DIMENSIONS: 335*510*180 mm

3.2 Installation

With VIVASALT Comfortable device is also added a paper template that allows you to accurately grab the position of the device on the wall and accurately measure the position of the holes into which the plastic screw inserts are inserted to secure the device to the wall (Figure 1).

The deflector shown in Figure 2 is used for mounting on the inner wall of the salt room.

The device is mounted on a wall where no other equipment is located within a distance of at least 50 cm from the device.



The device is mounted at a height of at least 85cm from the ground.

It is always mounted in a vertical (vertical) position. Usually, the device is located on the outer wall of the room into which the aerosol salt is blown in through a hole in the wall. When installing it, make sure that the installation is not visible, so that the device is not easily accessible to children or other unauthorized persons or visitors.

The device must be mounted in a way, that is reachable for personel and that the case can be opened.

Figure 1:

Marking the holes for the installation of the device using a template: Place the template on the wall and wipe the central aerosol blowing hole with a pencil and mark the holes through which the device is fastened with screws.



4. Working with VIVASALT Mediteran

4.1 Salt

The salt of high quality, which is pure, without additives and additives, is always used for the implementation of salt therapies. It is usually packaged in 25 kg bag and because it is free from additives, it wants to quickly harden. The salt is small (size 1 mm) and before being inserted into VIVASALT Mediteran it is necessary to make sure that there are no lumps, that they do not match and that the temperature is appropriate. Usually we break the grub with a kitchen mortar.

Salt Preparation Equipment:



KITCHEN MORTAR



INDUCTION PLATE




WOODEN KITCHEN SPOON



METAL POT (without cover)

The device may be damaged by inserting a salt that does not have a suitable temperature and if it contains clumps. If the salt is not dry enough, a salt crust is formed and the aerosol cant be made because the salt doesn't get to the mill.

A maximum of 30 g of salt may be added to the mill at a time (2 cups).





All information on how to handle the salt is obtained at the same time as training, which is usually done when the VIVASALT device is mounted.

Before starting the therapy, the parameters for performing the therapy should be set, that is, the duration of the program and the choice of the program.

The device must be connected to the power supply via a connecting cable. Connect the other end of the adapter cable to the AC adapter connector on the salt generator.



Before you press the "START" button, the parameters for the implementation of the therapy should be set.

The power cord must be plugged into the outlet. Switch the device on by pressing the switch located on the bottom of the device immediately next to the power cord and the Home screen is displayed (figure 1).



Figure 1: welcome screen



The screen is touch sensitive. Use a pen with a plastic tip to touch the screen and adjust the session parameters. Touching the green input arrow will open a new screen where the program is selected. Each of the 10 programs can be stored so that there is no need for additional storage and adjustment of the device during operation.

*** Save programs**

Pressing any of the suggested programs 1 to 10 will open a new image screen⁴. where you can adjust the selected program. Two status bars can be seen on the screen, one is for setting the session duration while the other bar can be used to adjust the power of the device. By touching the white fields of the status bars with the pen, you can adjust the time from 1-99 min and the power of the device from 1-100%.

When the values on both status bars are set, touching the “back” button in the upper left corner / green arrow return to the previous screen / will save the set parameters within the

selected program. Parameter storage can now be accessed in the same way on other device programs.

When all device programs have been saved, the device can be started.

5. Cleaning and maintenance of VIVASALT Mediteran

Maintenance and cleaning of the device is only carried out when it is disconnected from the main power. Never use wet cleaning or abrasive cleaning. The device should be completely cleaned at the end of the day. During the working day, after the end of each therapy, the dispenser is always emptied, all the salt that has remained, returns to the salt drying pot and cleans the mill with a brush.

5.1 Cleaning procedure

1. Turn off the device and unplug the power cord from the outlet
2. Unscrew the nut and gently pull the power cable connectors of the mill and the fan (figure 7, figure 7.1, figure 7.2)



Figure 7: Unscrew the nut

3. Carefully pull the elastic holders of the mill to the bottom and the side (Fig. 8). Pull them only as long as necessary so that the mill can be removed. If you get rough, you can cause deformation and may cause the device to fail. Carefully remove the mill in such a way that the cables of the device do not twist, tighten or break.




Figure 8: Removing the mill (magnets)

4. Carefully clean the mill using a brush (Fig. 9); the salt left in the mill belongs to the trash. Place it on the side, in a safe place until you are re-mounting it to the device.



Figure 9: Cleaning the mill

5. Remove the fan housing including the fan and dome by pulling



Figure 10: Removing the fan housing with the fan (magnets) and clean it with brush




Figure 11: Removing the dome and clean it with brush

6. Clean gently with the soft brush all the accessible areas of the fan and place the cleaned in a place that will not hinder the further cleaning of the device. Then re-mount it in its original position.
7. Clean all ventilation channels with a soft brush.
8. Check that the machine is well cleaned.
9. Return the cleaned parts of the device in the following order:

- a dosing device with a funnel
- fan housing with fan
- mill
- connect the connectors for the fan and mill
- Close the door.

The device is now ready for the next day.



6. Solving Problems

Error description	Probable cause of error	Solution
The device does not work	Check the fuses.	Replace the main fuse, check the connection of the cables, check on the bottom of the device if the fuse buttons of fan, mill and doser are convex and push them back.
	Problems with electrical voltage.	Check the power and mains cable.
The engine of the salt mill is loud or too hot	The salt mill is overloaded with salt or the engine is closed.	Clean the mill and make sure the salt meets the requirements (without additives and additives).
The salt mill works normally, but the salt does not get to it	Problems in the operation of the fan in the halogenerator and the quality of the salt.	Check the fan operation and check that the salt is not moist (it is best to dry it again).
Settings and programs do not appear when the halogener is turned on	The microprocessor memory is worn out	Reset programs or new upload is needed.
No aerosol.	Check the salt.	Disconnect the power cord from the outlet. Check if the salt is in the dispenser. Check if the salt is dry. Check if there is a failure of the fuse
Salt from the doser does not reach the salt mill	It is likely that there are larger lumps, a solid or moist salt, which prevents salt from entering the snail. The sanil is jammed, can not rotate, or the engine is dirty.	It is necessary to empty the dispenser and clean the snail. Dry the salt and crush the lumps (with the help of a mortar). Use only the recommended medical salt. Also, check that the connection are well connected.
Fuse failure on the device	Poorly cleaned parts of the mill, doser or fan.	Clean the machine exactly as described in the chapter "Cleaning", then press the fuse button located at the bottom of the device and start the therapy.



7. Warranty terms

VIVASALT devices are warranted if there is an error in the production itself or a material error:

VIVA SALT d.o.o. recognizes repair during the warranty period for unused products purchased from the manufacturer; the warranty is valid for 12 months from the date of purchase and is acknowledged for cases of materiel defects or manufacturing defects.

The guarantee is not recognized in the following cases:

1. The performance of the warranty is not possible in the following cases:
 - Damage to the device due to improper handling;
 - Improper maintenance and servicing
 - The halogenerator and the electrical installation were installed by an unauthorized person
 - If parts that are worn out or damaged are not replaced timely
 - In the event of irregular, incorrect cleaning or if the methods used are not recommended by the manufacturer
 - If an unauthorized person intervenes or is serviced by an unauthorized person.
2. The warranty shall be taken into account for errors arising from production or hidden material defects. Viva salt d.o.o. or his authorized representative will replace defective parts without additional costs for the buyer.
3. Viva salt d.o.o. or his representative does not accept liability for loss, damage or injury to persons or property, nor is there a guarantee of liability for cancellation or failure of the halogenerator and its installations if the work was carried out by an unauthorized person and without the approval of Viva salt d.o.o. or its representatives.
4. The warranty does not cover any costs or damages (business or non-business) that would arise in the event of failure of the device due to manufacturing defects or material.
5. The warranty does not apply to damage caused by a product not supplied by Viva salt d.o.o.
6. The warranty does not apply in case of damage caused by force majeure such as earthquake, floods, fire or lightning.
7. The warranty does not apply to damage caused by the supply voltage and its effects and irregularities in the existing installation.



8. The warranty is enforced in accordance with the following procedure:

- The buyer informs the authorized representative. He attempts to solve the problem via e-mail, phone, skype or other communication devices.

In case the problem is not resolved in this way:

- The buyer transmits, as agreed with the authorized representative, the defective part of the device for inspection and repair.
- The authorized representative corrects and returns to the customer the received part of the device. The warranty does not include service and maintenance of the device.
- The repair and transport are free of charge if it was found that the manufacturer's fault was.

7.1 Contraindications (responsibility to patients)

The buyer of the Viva salt device is responsible for its clients and will not mislead them and promise false results of treatment, especially if not recommended for people with the following health conditions:

- * severe hypertension (illness high blood pressure)
- * bleeding
- * Blood diseases in the acute phase
- * acute phase of respiratory diseases
- * severe chronic respiratory failure
- * active tuberculosis
- * mental illness and any consumption of illicit substances
- * cancerous diseases
- * poisoning
- * metabolic diseases
- * pregnancy
- * organ failure (especially liver and heart)
- * spitting blood
- * heart failure
- * chronic obstructive pulmonary disease

Salt therapy is not a substitute for medication and treatment. If the patient has a health problem, consult with your chosen doctor about the visit of salt therapies.

Important:

- **In case of doubt about the effect of salt therapies on health, the patient should consult with his doctor.**

8. Storage of the device

If it is necessary for any reason to disassemble the device, it must be stored in the original packaging in which the device was received. Before storing, the device needs to be cleaned well. Keep at temperatures between 5-30 degrees Celsius and humidity up to 60%.

